# AGPTEK

## FS10C

# User Manual

## Pulse Oximeter


# Content

**01/04  
05/08  
09/12  
13/16  
17/20**

**English  
Deutsch  
Italiano  
Español  
français**

# Introduction:


Pulse Oximeters are considered a suitable alternative, though by no means equal to, blood gas testing. It can be applied for Fitness Monitoring and High Altitude Sports. If you want to use for disease surveillance such as COPD Management, Asthma

Patient Care and Piece of Mind  
for Pulmonary Patients, please  
consult your doctor.

## How to Use a Pulse Oximeter Properly

1. Press the bottom to open the probe.

2. Insert the  
finger.

3. Press button to activate the  
Oximeter, you will need to wait  
for 5 more seconds for the Pulse

Oximeter to stabilize its readings.  
**To get the best reading possible,  
you should NOT:**

- Test with your painted fingers and toenails which can obstruct proper readings – use an earlobe instead in this case.
- Test under bright overhead lighting, including sunlight can effect readings – take your readings in rooms with normal or low-light conditions.
- Move while testing, body movement including shivering, sneezing, etc., can all affect proper readings.

**Normal Range for Healthy  
Individuals in a Quiet State:  
SpO2% >94      60<PR<120**

For detailed instructions in  
English, please visit  
<http://www.agptek.com>.

If you have any problem or  
suggestions,please contact  
[support@agptek.com](mailto:support@agptek.com).

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