AGPTEK

User Manual

Model: AS-302-L

EN > DE > FR > IT > ES

English

01 - 05

Deutsch

06 - 11

Français

12 - 17

Italiano


18 - 23

Español

24 - 29


Introduction

Pulse Oximeters are considered a suitable alternative, though by no means equal to, blood gas testing. It can be applied for Fitness Monitoring and High Altitude Sports. If you want to use for disease surveillance such as COPD Management, Asthma Patient Care and Piece of Mind for Pulmonary Patients, please consult your doctor.

- ① Low Power Indication
- ② SpO2 Value
- ③ Pulse Rate Value
- ④ Pulse Bar Graph
- ⑤ Operating Button

How to Use a Pulse Oximeter Properly


1. Press the bottom to open the probe.
2. Insert the finger.
3. Press button **5** to activate the Oximeter, you will need to wait for 8 more seconds for the Pulse Oximeter to stabilize its readings.

To get the best reading possible, you should NOT:

- Test with your painted fingers and toenails which can obstruct proper readings - use an earlobe instead in this case.
- Test under bright overhead lighting including sunlight, can effect readings - take your readings in rooms with normal or low-light conditions.
- Move while testing, body movement including shivering, sneezing, etc.,

can all affect proper readings.

Battery Cover Installation

Normal Range for Healthy Individuals in a Quiet State:

SpO₂% >94

60 < PR < 120

If you have any problem or suggestion, please contact support@agptek.com.

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