AGPTEK

Back Posture Corrector

# CONTENT

English	01-04
Deutsch	05-08
Français	09-12
Italiano	13-16
Español	17-20
日本語	21-24

Thank you for purchasing AGPTEK Back Posture Corrector. Our product is mainly used to correct your posture by stretching your shoulders and fixing your shoulder blades.

## > To our customers

ⓘ **Attention!** It is recommended that the user should be over 12 years old and wear it for no less than 30 minutes every time. If you feel terrible, please take it off and relax yourself.

Back posture corrector

Multifunctional pad

## > How to wear


1. Take posture corrector out. Adjust its belts to the bottom of shoulder straps.

2. Slip your arms through each side of the posture brace.

3. Put on the corrector as if you wear a vest.


4. Adjust hook & loop fasteners and front tightening straps according to your needs.

## ➤ Instructions for loosened straps

### ☆ Without a partner

☆ With a partner

If the straps is loosened, please refer to the following steps to wear it with the help of a partner.

1. Place the straps on your shoulder, two belts falling on the front of your belly.

2. Your partner helps you put the straps through buckles and pull them out.

3. Take belts from your partner' s hand.

4. Pull on straps slowly and evenly to correct posture.

5. Done!

## ➤ How to use the multifunctional pad.

1. After wearing the corrector, place the multifunctional pad on your belly.

2. Attach one belt to the middle of pad.

3. The other belt also attaches to the middle of pad.

4. Roll up the pad in a bottom-up direction.

5. Fold the top part of pad (with hook and loop fasteners) to the pad that has been rolled. Make the belts wrapped by the pad tightly.

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