AGPTEK

FS20D

User Manual


Pulse Oximeter

Content

English	01/04
Deutsch	05/09
Italiano	10/14
Español	15/18

Introduction:


Pulse Oximeters are considered a suitable alternative, though by no means equal to, blood gas testing. It can be applied for Fitness Monitoring and High Altitude Sports. If you want to use for disease surveillance such as COPD Management, Asthma


Patient Care and Piece of Mind for Pulmonary Patients, please consult your doctor.

How to Use a Pulse Oximeter Properly

1. Press the bottom to open the probe.

2. Insert the finger.

3. Press button  to activate the Oximeter, you will need to wait for 5 more seconds for the Pulse

Oximeter to stabilize its readings.
To get the best reading possible, you should NOT:

- Test with your painted fingers and toenails which can obstruct proper readings – use an earlobe instead in this case.
- Test under bright overhead lighting, including sunlight can affect readings – take your readings in rooms with normal or low-light conditions.
- Move while testing, body movement including shivering, sneezing, etc., can all affect proper readings.

Normal Range for Healthy Individuals in a Quiet State:

SpO₂% >94 60 < PR < 120

When the measurement is beyond the normal range (SpO₂% < 90 or 40 < PR < 120), the Oximeter will alarm which probably means that your body may have a problem.

For detailed instructions in English, please visit <http://www.agptek.com>.

If you have any problem or suggestions, please contact support@agptek.com.

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