
AGPTEK

Welcome Guide

Model: YM101

EN > DE > FR > IT > ES

English

01 - 06

Deutsch

07 - 13

Français

14 - 20

Italiano


21 - 27

Español

28 - 34


Introduction


Pulse Oximeters are considered a suitable alternative, though by no means equal to, blood gas testing. It can be applied for Fitness Monitoring and High Altitude Sports. If you want to use for disease surveillance such as COPD Management, Asthma Patient Care and Piece of Mind for Pulmonary Patients, please consult your doctor.

- ① SpO2 Value
- ② Operating Button
- ③ Low Power Indication
- ④ Pulse Rate Value
- ⑤ Pulse Bar Graph
- ⑥ Perfusion Index

How to Use a Pulse Oximeter Properly

1. Install batteries.
2. Press the bottom to open the probe.
3. Insert the finger.
4. Press button **2** to activate the oximeter, you will need to wait for 8 more seconds for the pulse oximeter to stabilize its readings. (Note: Support showing readings in

2 directions, short press to change the reading direction.)

Customize Settings

Non-professionals are advised not to customize the settings and use the default mode.

After power on, long press button for entering Settings, short press to view the value of each item, and long press to enter the item's submenu. Then, short press to select the parameter and long press to confirm.

Item 1: Setup the brightness of LED display

Item 2: Setup the SpO2 alarm limits

Item 3: Setup the PR alarm limits

Item 4: Turn alarm on/off

Item 5: Check the software version

To get the best reading possible, you should NOT:

- Test with your painted fingers and toenails which can obstruct proper readings - use an earlobe instead in this case.
- Test under bright overhead lighting including sunlight, can effect readings - take your readings in rooms with normal or low-light conditions.
- Move while testing, body

movement including shivering, sneezing, etc., can all affect proper readings.

Normal Range for Healthy Individuals in a Quiet State:

SpO₂% >94

60 < PR < 120

For more instructions, please see the detailed user manual.

If you have any problem or suggestion, please contact support@agptek.com.

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